



TELO2 rabbit pAb

Catalog No	YP-Ab-09892
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TELO2 KIAA0683
Immunogen	Synthesized peptide derived from human TELO2 AA range: 605-655
Specificity	This antibody detects endogenous levels of human TELO2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TEL2;Protein clk-2 homolog;Protein clk 2 homolog;KIAA0683;hCLK2;TELO2;TELO2 KIAA0683;Telomere length regulation protein TEL2 homolog (Protein clk-2 homolog) (hCLK2)
Observed Band	92 kDa
Calculated Molecular Weight	837 aa, 92 kDa
Cell Pathway	Cytoplasm. Membrane. Nucleus. Chromosome, telomere .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
