



GPI8 rabbit pAb

Catalog No	YP-Ab-09902
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PIGK GPI8
Immunogen	Synthesized peptide derived from human GPI8 AA range: 331-381
Specificity	This antibody detects endogenous levels of human GPI8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Phosphatidylinositol-glycan biosynthesis class K protein;hGPI8;GPI-anchor transamidase;GPI8;GPI transamidase;PIGK GPI8;GPI-anchor transamidase (GPI transamidase) (EC 3.-.-.-) (GPI8 homolog) (hGPI8) (Phosphatidylinositol-glycan biosynthesis class K protein) (PIG-K)
Observed Band	40-45 kDa
Calculated Molecular Weight	45 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

