



LEMD1 rabbit pAb

Catalog No	YP-Ab-09906
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LEMD1
Immunogen	Synthesized peptide derived from human LEMD1 AA range: 12-62
Specificity	This antibody detects endogenous levels of human LEMD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LEMD1;LEM domain-containing protein 1 (Cancer/testis antigen 50) (CT50) (LEM domain protein 1) (LEMP-1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Testis-specific. Isoform 6 is detected in 17 of 18 colon cancer tissues examined.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
