



VEZA rabbit pAb

Catalog No	YP-Ab-09927
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VEZT
Immunogen	Synthesized peptide derived from human VEZA AA range: 325-375
Specificity	This antibody detects endogenous levels of human VEZA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VEZA;VEZT;Vezatin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cell projection, stereocilium membrane . Cell junction, adherens junction . Nucleus . Cytoplasmic vesicle, secretory vesicle, acrosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
