



SKAP rabbit pAb

Catalog No	YP-Ab-09930
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SKAP C15orf23 TRAF4AF1 HSD11
Immunogen	Synthesized peptide derived from human SKAP AA range: 34-84
Specificity	This antibody detects endogenous levels of human SKAP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SKAP;SKAP C15orf23 TRAF4AF1 HSD11;Small kinetochore-associated protein (SKAP) (Kinetochore-localized astrin-binding protein) (Kinastrin) (TRAF4-associated factor 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome, centromere, kinetochore . Cytoplasm, cytoskeleton, spindle pole . Cytoplasm, cytoskeleton, microtubule organizing center . Colocalizes with microtubules around centrosomes in prophase and with the mitotic spindle at prometaphase and metaphase. From late prometaphase to anaphase, is highly concentrated on kinetochores. Located at the kinetochore-microtubule interface. The astrin (SPAG5)-kinastrin (SKAP) complex localizes to the microtubule plus ends (PubMed:23035123). .
Tissue Specificity	Widely expressed, including in skin.
Function	
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images