



bRHBT1 rabbit pA

Catalog No	YP-Ab-09933
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RHOBTB1 KIAA0740
Immunogen	Synthesized peptide derived from human RHBT1 AA range: 296-346
Specificity	This antibody detects endogenous levels of human RHBT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RHBT1;RHOBTB1 KIAA0740;Rho-related BTB domain-containing protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Ubiquitous, with highest levels in skeletal muscle, placenta, testis, stomach, and kidney, followed by uterus and adrenal gland. Expressed in a variety of fetal tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
