



MF2L2 rabbit pAb

Catalog No	YP-Ab-09934
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MCF2L2 DRG KIAA0861
Immunogen	Synthesized peptide derived from human MF2L2 AA range: 97-147
Specificity	This antibody detects endogenous levels of human MF2L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MF2L2;MCF2L2 DRG KIAA0861;Probable guanine nucleotide exchange factor MCF2L2 (Dbs-related Rho family guanine nucleotide exchange factor) (MCF2-transforming sequence-like protein 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Significantly expressed in brain and modestly in pancreas, brain and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

