



CSMD1 rabbit pAb

Catalog No	YP-Ab-09936
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CSMD1 KIAA1890 UNQ5952/PRO19863
Immunogen	Synthesized peptide derived from human CSMD1 AA range: 5-55
Specificity	This antibody detects endogenous levels of human CSMD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CSMD1;CSMD1 KIAA1890 UNQ5952/PRO19863;CUB and sushi domain-containing protein 1 (CUB and sushi multiple domains protein 1)
Observed Band	392kD
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Weakly expressed in most tissues, except in brain. Expressed at intermediate level in brain, including cerebellum, substantia nigra, hippocampus and fetal brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

