



SN12L rabbit pAb

Catalog No	YP-Ab-09941
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SLFN12L SLFN5
Immunogen	Synthesized peptide derived from human SN12L AA range: 147-197
Specificity	This antibody detects endogenous levels of human SN12L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SN12L;SLFN12L SLFN5;Schlafen family member 12-like
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane; Single-pass membrane protein.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
