



RN19B rabbit pAb

Catalog No	YP-Ab-09943
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RNF19B IBRDC3 NKLAM
Immunogen	Synthesized peptide derived from human RN19B AA range: 272-322
Specificity	This antibody detects endogenous levels of human RN19B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RN19B;RNF19B IBRDC3 NKLAM;E3 ubiquitin-protein ligase RNF19B (EC 6.3.2.-) (IBR domain-containing protein 3) (Natural killer lytic-associated molecule) (RING finger protein 19B)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic granule membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed specifically in natural killer cells, activated macrophages and cytotoxic T-cells (PubMed:10438909). Present in natural killer cells (at protein level) (PubMed:10912506). Ubiquitously expressed with high expression in testis (PubMed:27485036).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images