



NBPFF rabbit pAb

Catalog No	YP-Ab-09947
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	NBPF15
Immunogen	Synthesized peptide derived from human NBPFF AA range: 382-432
Specificity	This antibody detects endogenous levels of human NBPFF
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NBPFF;NBPF15;Neuroblastoma breakpoint family member 15
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitously expressed with a higher expression observed in breast and liver. Also expressed in neuroblastoma cell line.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

