



GPDM rabbit pAb

Catalog No	YP-Ab-09954
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPD2
Immunogen	Synthesized peptide derived from human GPDM AA range: 310-360
Specificity	This antibody detects endogenous levels of human GPDM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.1.5.3;GDH2;Glycerol-3-phosphate dehydrogenase;mitochondrial;GPD 2; GPD M;GPDM;GPD2;mitochondrial (GPD-M) (GPDH-M) (EC 1.1.5.3) (mtGPD)
Observed Band	68 kDa
Calculated Molecular Weight	41 kDa, 81 kDa
Cell Pathway	Mitochondrion.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

