



BI1 rabbit pAb

Catalog No	YP-Ab-09958
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMBIM6 BI1 TEGT
Immunogen	Synthesized peptide derived from human BI1 AA range: 93-143
Specificity	This antibody detects endogenous levels of human BI1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Bax inhibitor 1;BAXI1;BI 1;BI1;TEGT;TMBIM6;TMBIM6 BI1 TEGT;Bax inhibitor 1 (BI-1) (Testis-enhanced gene transcript protein) (Transmembrane BAX inhibitor motif-containing protein 6)
Observed Band	27 kDa
Calculated Molecular Weight	27 kDa
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Highly abundant in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
