



YRDC rabbit pAb

Catalog No	YP-Ab-09962
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	YRDC DRIP3 IRIP
Immunogen	Synthesized peptide derived from human YRDC AA range: 214-264
Specificity	This antibody detects endogenous levels of human YRDC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DRIP3;hIRIP;IRIP;RP11 109P14.4;SUA5;YRDC;YRDC DRIP3 IRIP;YrdC domain-containing protein;mitochondrial (Dopamine receptor-interacting protein 3) (Ischemia/reperfusion-inducible protein homolog) (hIRIP)
Observed Band	29 kDa
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Mitochondrion . Cell membrane ; Peripheral membrane protein . A large fraction localizes in the cytoplasm, whereas a smaller fraction is imported to mitochondria. .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
