



PAG15 rabbit pAb

Catalog No	YP-Ab-09972
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLA2G15 LYPLA3 UNQ341/PRO540
Immunogen	Synthesized peptide derived from human PAG15 AA range: 233-283
Specificity	This antibody detects endogenous levels of human PAG15
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAG15;PLA2G15 LYPLA3 UNQ341/PRO540;Group XV phospholipase A2 (EC 2.3.1.-) (1-O-acylceramide synthase) (ACS) (LCAT-like lysophospholipase) (LLPL) (Lysophospholipase 3) (Lysosomal phospholipase A2) (LPLA2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Lysosome . Secreted . Membrane ; Peripheral membrane protein .
Tissue Specificity	Detected in blood plasma (at protein level) (PubMed:10092508, PubMed:20410020). Ubiquitous. Highly expressed in heart, placenta, skeletal muscle, kidney and pancreas. Detected at lower levels in spleen, thymus, prostate, testis, ovary, small intestine, colon and peripheral blood leukocytes (PubMed:10092508).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images