



DUS8 rabbit pAb

Catalog No	YP-Ab-09999
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DUSP8 C11orf81 VH5
Immunogen	Synthesized peptide derived from human DUS8 AA range: 98-148
Specificity	This antibody detects endogenous levels of human DUS8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C11orf81;Chromosome 11 open reading frame 81;Deubiquitinating enzyme 8; Dual Specificity Phosphatase 8;Dual specificity phosphatase M3/6;Dual specificity protein phosphatase 8;Dual specificity protein phosphatase hVH 5; Dual specificity protein phosphatase hVH-5;DUS8;DUS8_HUMAN;DUSP8;EC 3.1.3.16;EC 3.1.3.48;FLJ42476;FLJ42958;H1 phosphatase vaccinia virus homolog;H1 phosphatase;vaccinia virus;homolog of;HB5;hUBPy;HVVH5;HVVH8; M3/6;mouse;Neuronal tyrosine threonine phosphatase 1;Serine/threonine specific protein phosphatase;Ubiquitin carboxyl terminal hydrolase 8;Ubiquitin specific processing protease 8;Ubiquitin thiolesterase 8;VH5;DUSP8 C11orf81 VH5;Dual specificity protein phosphatase 8 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH-5)
Calculated Molecular Weight	69kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Abundant in brain, heart and skeletal muscle.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

