



ZDHC1 rabbit pAb

Catalog No	YP-Ab-10011
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ZDHC1 C16orf1 ZNF377
Immunogen	Synthesized peptide derived from human ZDHC1 AA range: 81-131
Specificity	This antibody detects endogenous levels of human ZDHC1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZDHC1;ZDHC1 C16orf1 ZNF377;Probable palmitoyltransferase ZDHC1 (EC 2.3.1.-) (DHHC domain-containing cysteine-rich protein 1) (Zinc finger DHHC domain-containing protein 1) (DHHC-1) (Zinc finger protein 377)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endosome membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus .
Tissue Specificity	Widely expressed with significant expression in heart, brain, placenta, lung, liver, kidney, testis, thymus and small intestine (PubMed:16647879). Expressed at lower levels in adult pancreas and lung (PubMed:10395086).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images