



VOPP1 rabbit pAb

Catalog No	YP-Ab-10017
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VOPP1 ECOP GASP
Immunogen	Synthesized peptide derived from human VOPP1 AA range: 120-170
Specificity	This antibody detects endogenous levels of human VOPP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ECOP;GASP;VOPP1;VOPP1 ECOP GASP;Vesicular;overexpressed in cancer; prosurvival protein 1 (EGFR-coamplified and overexpressed protein) (ECop) (Glioblastoma-amplified secreted protein) (Putative NF-kappa-B-activating protein 055N)
Observed Band	
Calculated Molecular Weight	172 aa, 19 kDa
Cell Pathway	Cytoplasmic vesicle membrane ; Single-pass type I membrane protein . When overexpressed, localizes in the nucleus and perinuclear regions.
Tissue Specificity	Widely expressed with highest levels in thymus and ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
