



# CBP 35 Monoclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                                                                                                                              |
|------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-00350                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Applications</b>                | WB                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Gene Name</b>                   | LGALS3                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human Galectin 3. AA range:141-190                                                                                                                                                                                                                                                                                                                       |
| <b>Specificity</b>                 | CBP 35 Monoclonal Antibody detects endogenous levels of CBP 35 protein.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                                                         |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Synonyms</b>                    | 35 kDa lectin;Carbohydrate-binding protein 35;CBP 35;CBP35;Gal 3;Galectin-3;LGALS3;MAC2;Gal-3;Galactose-specific lectin 3;Galactoside-binding protein;GALBP;IgE-binding protein;L-31;Laminin-binding protein;Lectin L-29;Mac-2 antigen                                                                                                                                                                                       |
| <b>Observed Band</b>               | 33kD                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Calculated Molecular Weight</b> | 33kD                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Cell Pathway</b>                | Cytoplasm . Nucleus. Secreted . Secreted by a non-classical secretory pathway and associates with the cell surface. Can be secreted; the secretion is dependent on protein unfolding and facilitated by the cargo receptor TMED10; it results in protein translocation from the cytoplasm into the ERGIC (endoplasmic reticulum-Golgi intermediate compartment) followed by vesicle entry and secretion (PubMed:32272059). . |
| <b>Tissue Specificity</b>          | A major expression is found in the colonic epithelium. It is also abundant in the activated macrophages. Expressed in fetal membranes.                                                                                                                                                                                                                                                                                       |
| <b>Function</b>                    | function:Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimulation by CSPG4 of endothelial cells migration. Together with DMBT1, required for terminal differentiation of columnar epithelial cells during early embryogenesis.,online information:Galectin-3,similarity:Contains 1                                                                                           |



galectin domain.,subcellular location:Cytoplasmic in adenomas and carcinomas. May be secreted by a non-classical secretory pathway and associate with the cell surface.,subunit:Probably forms homo- or heterodimers. Interacts with DMBT1 (By similarity). Forms a complex with the ITGA3, ITGB1 and CSPG4. Interacts with LGALS3BP, LYPD3, CYHR1 and UACA.,tissue specificity:A major expression is found in the colonic epithelium. It is also abundant in the activated macrophages.,

#### Background

This gene encodes a member of the galectin family of carbohydrate binding proteins. Members of this protein family have an affinity for beta-galactosides. The encoded protein is characterized by an N-terminal proline-rich tandem repeat domain and a single C-terminal carbohydrate recognition domain. This protein can self-associate through the N-terminal domain allowing it to bind to multivalent saccharide ligands. This protein localizes to the extracellular matrix, the cytoplasm and the nucleus. This protein plays a role in numerous cellular functions including apoptosis, innate immunity, cell adhesion and T-cell regulation. The protein exhibits antimicrobial activity against bacteria and fungi. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Oct 2014],


#### matters needing attention

Avoid repeated freezing and thawing!

#### Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images



Western Blot analysis of various cells using CBP 35 Monoclonal Antibody