



CSRNP3 Monoclonal Antibody

Catalog No	YP-mAb-00362
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CSRNP3
Immunogen	The antiserum was produced against synthesized peptide derived from human TAIP-2. AA range:171-220
Specificity	CSRNP3 Monoclonal Antibody detects endogenous levels of CSRNP3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CSRNP3;Cysteine/serine-rich nuclear protein 3;FAM130A2;TAIP2;CSRNP-3;Protein FAM130A2;TGF-beta-induced apoptosis protein 2;TAIP-2
Observed Band	64kD
Calculated Molecular Weight	64kD
Cell Pathway	Nucleus .
Tissue Specificity	Embryo,Thymus,
Function	function:Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis.,sequence caution:Unusual initiator. The initiator methionine is coded by a non-canonical CTG leucine codon.,similarity:Belongs to the AXUD1 family.,
Background	function:Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis.,sequence caution:Unusual initiator. The initiator methionine is coded by a non-canonical CTG leucine codon.,similarity:Belongs to the AXUD1 family.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CSRNP3 Monoclonal Antibody