



DAP-5 Monoclonal Antibody

Catalog No	YP-mAb-00370
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EIF4G2
Immunogen	The antiserum was produced against synthesized peptide derived from human EIF4G2. AA range:41-90
Specificity	DAP-5 Monoclonal Antibody detects endogenous levels of DAP-5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DAP-5;Disks large homolog 7;Disks large-associated protein 5;Disks large-associated protein DLG7;EIF4G2;Eukaryotic translation initiation factor 4 gamma 2;DAP5;OK/SW-cl.75;eIF-4-gamma 2;eIF-4G 2;eIF4G 2;Death-associated protein 5;p97
Observed Band	90kD
Calculated Molecular Weight	90kD
Cell Pathway	cytosol,cell-cell adherens junction,membrane,eukaryotic translation initiation factor 4F complex,axon,
Tissue Specificity	Ubiquitously expressed in all adult tissues examined, with high levels in skeletal muscle and heart. Also expressed in fetal brain, lung, liver and kidney.
Function	function:Appears to play a role in the switch from cap-dependent to IRES-mediated translation during mitosis, apoptosis and viral infection. Cleaved by some caspases and viral proteases.,miscellaneous:This gene has been shown to be extensively edited in the liver of APOBEC1 transgenic animal model. Its aberrant editing could contribute to the potent oncogenesis induced by overexpression of APOBEC1. The aberrant edited sequence, called NAT1, is likely to be a fundamental translational repressor.,PTM:Phosphorylation; hyperphosphorylated during mitosis.,similarity:Belongs to the eIF4G family.,



similarity:Contains 1 MI domain.,similarity:Contains 1 MIF4G domain.,
similarity:Contains 1 W2 domain.,subunit:Interacts with the serine/threonine
protein kinases MKNK1 and MKNK2. Binds EIF4A and EIF3. Interacts with
MIF4GD.,tissue specificity:Ubiquitously expressed in all adult tissues examined,
with h

Background

Translation initiation is mediated by specific recognition of the cap structure by eukaryotic translation initiation factor 4F (eIF4F), which is a cap binding protein complex that consists of three subunits: eIF4A, eIF4E and eIF4G. The protein encoded by this gene shares similarity with the C-terminal region of eIF4G that contains the binding sites for eIF4A and eIF3; eIF4G, in addition, contains a binding site for eIF4E at the N-terminus. Unlike eIF4G, which supports cap-dependent and independent translation, this gene product functions as a general repressor of translation by forming translationally inactive complexes. In vitro and in vivo studies indicate that translation of this mRNA initiates exclusively at a non-AUG (GUG) codon. Alternatively spliced transcript variants encoding different isoforms of this gene have been described. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using DAP-5 Monoclonal Antibody