



FASTKD1 Monoclonal Antibody


Catalog No	YP-mAb-00396
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FASTKD1
Immunogen	The antiserum was produced against synthesized peptide derived from human FASTKD1. AA range:561-610
Specificity	FASTKD1 Monoclonal Antibody detects endogenous levels of FASTKD1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FASTKD1;FAST kinase domain-containing protein 1;KIAA1800
Observed Band	97kD
Calculated Molecular Weight	97kD
Cell Pathway	Mitochondrion . Preferentially localizes to mitochondrial RNA granules, platforms for post-transcriptional RNA modification and ribosome assembly (PubMed:28335001). .
Tissue Specificity	Expression detected in spleen, thymus, testis, ovary, colon, heart, smooth muscle, kidney, brain, lung, liver and white adipose tissue with highest expression in heart.
Function	caution:It is uncertain whether Met-1 or Met-15 is the initiator.,similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,
Background	caution:It is uncertain whether Met-1 or Met-15 is the initiator.,similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using FASTKD1 Monoclonal Antibody