



FASTKD3 Monoclonal Antibody

Catalog No	YP-mAb-00398
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FASTKD3
Immunogen	The antiserum was produced against synthesized peptide derived from human FASTKD3. AA range:121-170
Specificity	FASTKD3 Monoclonal Antibody detects endogenous levels of FASTKD3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FASTKD3;FAST kinase domain-containing protein 3
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Expression detected in spleen, thymus, testis, ovary, colon, heart, smooth muscle, kidney, brain, lung, liver and white adipose tissue with highest expression in liver and thyroid.
Function	similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,
Background	This gene encodes a member of a small family of Fas-activated serine/threonine kinase domain (FASTKD) containing proteins that share an amino terminal mitochondrial targeting domain and multiple carboxy terminal FAST domains as well as a putative RNA-binding RAP domain. The members of this family are ubiquitously expressed and are generally most abundant in mitochondria-enriched tissues such as heart, skeletal muscle and brown-adipose tissue. Some members of this protein family may play a role in apoptosis. The protein encoded by this gene interacts with components of the mitochondrial respiratory



and translation networks. A pseudogene of this gene is also present on chromosome 5. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using FASTKD3 Monoclonal Antibody