



GNL3L Monoclonal Antibody

Catalog No	YP-mAb-00414
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GNL3L
Immunogen	The antiserum was produced against synthesized peptide derived from human GNL3L. AA range:71-120
Specificity	GNL3L Monoclonal Antibody detects endogenous levels of GNL3L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GNL3L;Guanine nucleotide-binding protein-like 3-like protein
Observed Band	66kD
Calculated Molecular Weight	66kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Lung,Skin,
Function	domain:In contrast to other GTP-binding proteins, this family is characterized by a circular permutation of the GTPase motifs described by a G4-G1-G3 pattern., function:Required for normal processing of ribosomal pre-rRNA. Required for cell proliferation. Binds GTP.,similarity:Belongs to the MMR1/HSR1 GTP-binding protein family.,
Background	The protein encoded by this gene appears to be a nucleolar GTPase that is essential for ribosomal pre-rRNA processing and cell proliferation. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, May 2010],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using GNL3L Monoclonal Antibody