




LDOC1L Monoclonal Antibody

Catalog No	YP-mAb-00432
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LDOC1L
Immunogen	The antiserum was produced against synthesized peptide derived from human LDOC1L. AA range:141-190
Specificity	LDOC1L Monoclonal Antibody detects endogenous levels of LDOC1L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LDOC1L;Protein LDOC1L;MAR6;MART6;Leucine zipper protein down-regulated in cancer cells-like;Mammalian retrotransposon-derived protein 6
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	
Tissue Specificity	Amygdala,Brain,Trachea,
Function	similarity:Belongs to the LDOC1 family.,
Background	similarity:Belongs to the LDOC1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LDOC1L Monoclonal Antibody