




NSE2 Monoclonal Antibody

Catalog No	YP-mAb-00465
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM84B
Immunogen	The antiserum was produced against synthesized peptide derived from human FAM84B. AA range:11-60
Specificity	NSE2 Monoclonal Antibody detects endogenous levels of NSE2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRAT domain-containing 2;LRATD2;NSE2;Protein FAM84B;Protein LRATD2; FAM84B;BCMP101;Breast cancer membrane protein 101;Protein NSE2
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	cytoplasm,plasma membrane,
Tissue Specificity	Expressed in esophageal squamous cell carcinomas.
Function	similarity:Belongs to the FAM84 family.,tissue specificity:Expressed in esophageal squamous cell carcinomas.,
Background	similarity:Belongs to the FAM84 family.,tissue specificity:Expressed in esophageal squamous cell carcinomas.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NSE2 Monoclonal Antibody