




ORAOV1 Monoclonal Antibody

Catalog No	YP-mAb-00469
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ORAOV1
Immunogen	The antiserum was produced against synthesized peptide derived from human ORAV1. AA range:81-130
Specificity	ORAOV1 Monoclonal Antibody detects endogenous levels of ORAOV1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CIAB1;LTO1;Oral cancer-overexpressed protein 1;Protein LTO1 homolog;TAOS1;ORAOV1;Tumor-amplified and overexpressed sequence 1
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed. Highly expressed in placenta, kidney and skeletal muscle.
Function	miscellaneous:Amplified and overexpressed in oral cancer cells., similarity:Belongs to the ORAOV1 family.,tissue specificity:Widely expressed. Highly expressed in placenta, kidney and skeletal muscle.,
Background	miscellaneous:Amplified and overexpressed in oral cancer cells., similarity:Belongs to the ORAOV1 family.,tissue specificity:Widely expressed. Highly expressed in placenta, kidney and skeletal muscle.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ORAOV1 Monoclonal Antibody