



PDRG1 Monoclonal Antibody

Catalog No	YP-mAb-00489
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PDRG1
Immunogen	The antiserum was produced against synthesized peptide derived from human PDRG1. AA range:58-107
Specificity	PDRG1 Monoclonal Antibody detects endogenous levels of PDRG1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C20orf126;p53 and DNA damage regulated 1;p53 and DNA damage-regulated protein 1;PDRG;PDRG1
Observed Band	19kD
Calculated Molecular Weight	19kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Predominantly expressed in normal testis and exhibits reduced but detectable expression in other organs.
Function	function:May play a role in chaperone-mediated protein folding .,induction:By UV and repressed by TP53/p53.,similarity:Belongs to the prefoldin subunit beta family.,tissue specificity:Predominantly expressed in normal testis and exhibits reduced but detectable expression in other organs.,
Background	function:May play a role in chaperone-mediated protein folding .,induction:By UV and repressed by TP53/p53.,similarity:Belongs to the prefoldin subunit beta family.,tissue specificity:Predominantly expressed in normal testis and exhibits reduced but detectable expression in other organs.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using PDRG1 Monoclonal Antibody