




Serinc1 Monoclonal Antibody

Catalog No	YP-mAb-00521
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SERINC1
Immunogen	The antiserum was produced against synthesized peptide derived from human SERC1. AA range:381-430
Specificity	Serinc1 Monoclonal Antibody detects endogenous levels of Serinc1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1253;SERINC1;serine incorporator 1;TDE1L;TDE2;TMS 2;Tumor differentially expressed protein 1-like;Tumor differentially expressed protein 2
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Adrenal gland,Brain,Epithelium,Liver,Testis,
Function	function:Enhances the incorporation of serine into phosphatidylserine and sphingolipids.,similarity:Belongs to the TDE1 family.,subunit:Interacts with SPTLC1.,
Background	function:Enhances the incorporation of serine into phosphatidylserine and sphingolipids.,similarity:Belongs to the TDE1 family.,subunit:Interacts with SPTLC1.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Serinc1 Monoclonal Antibody