




TDE2L Monoclonal Antibody

Catalog No	YP-mAb-00532
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SERINC2
Immunogen	The antiserum was produced against synthesized peptide derived from human SERC2. AA range:341-390
Specificity	TDE2L Monoclonal Antibody detects endogenous levels of TDE2L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAI;PAI2;PAI-2;PLANH2;HsT1201;SERPINB2;TDE2L;SERINC2;Serine incorporator 2;FKSG84;Tumor differentially expressed protein 2-like
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Colon,Liver,Pancreas,Uterus,
Function	similarity:Belongs to the TDE1 family.,
Background	similarity:Belongs to the TDE1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TDE2L Monoclonal Antibody