



TP53INP1 Monoclonal Antibody

Catalog No	YP-mAb-00535
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TP53INP1
Immunogen	The antiserum was produced against synthesized peptide derived from human TP53INP1. AA range:171-220
Specificity	TP53INP1 Monoclonal Antibody detects endogenous levels of TP53INP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	p53DINP1;SIP;Stress induced protein;Teap;TP53DINP1;TP53INP1;TP53INP1A;TP53INP1B;Tumor protein p53-inducible nuclear protein 1;Stress-induced protein;p53-dependent damage-inducible nuclear protein 1
Observed Band	27kD
Calculated Molecular Weight	27kD
Cell Pathway	Cytoplasm, cytosol. Nucleus. Nucleus, PML body. Cytoplasmic vesicle, autophagosome. Shuttles between the nucleus and the cytoplasm, depending on cellular stress conditions, and re-localizes to autophagosomes on autophagy activation.
Tissue Specificity	Ubiquitously expressed.
Function	function:In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis.,induction:By adriamycin, gamma irradiation and H ₂ O ₂ , in a p53-dependent way. At lower levels by ultraviolet irradiation.,subcellular location:Concentrated in PML/POD/ND10 nuclear bodies.,subunit:Interacts with TP53 and HIPK2.,tissue specificity:Ubiquitously expressed.,



Background

function: In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis., induction: By adriamycin, gamma irradiation and H₂O₂, in a p53-dependent way. At lower levels by ultraviolet irradiation., subcellular location: Concentrated in PML/POD/ND10 nuclear bodies., subunit: Interacts with TP53 and HIPK2., tissue specificity: Ubiquitously expressed.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using TP53INP1 Monoclonal Antibody