



ABCD4 Monoclonal Antibody

Catalog No	YP-mAb-00656
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABCD4
Immunogen	The antiserum was produced against synthesized peptide derived from human ABCD4. AA range:111-160
Specificity	ABCD4 Monoclonal Antibody detects endogenous levels of ABCD4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABC41;ABCD4;EST352188;P70R;P79R;PMP69;PMP70 related protein;PXMP1 L;PXMP1L;ATP-binding cassette sub-family D member 4;PMP70-related protein; Peroxisomal membrane protein 1-like;PXMP1-L;Peroxisomal membrane protein 69
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Lysosome membrane ; Multi-pass membrane protein . Targeted by LMBRD1 lysosomal chaperone to the lysosomal membrane. .
Tissue Specificity	Ubiquitous.
Function	similarity:Belongs to the ABC transporter family. ALD subfamily., similarity:Contains 1 ABC transmembrane type-1 domain.,similarity:Contains 1 ABC transporter domain.,subunit:Homodimer or heterodimer .,tissue specificity:Ubiquitous.,



Background

The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the ALD subfamily, which is involved in peroxisomal import of fatty acids and/or fatty acyl-CoAs in the organelle. All known peroxisomal ABC transporters are half transporters which require a partner half transporter molecule to form a functional homodimeric or heterodimeric transporter. The function of this peroxisomal membrane protein is unknown. However, it is speculated that it may function as a heterodimer for another peroxisomal ABC transporter and, therefore, may modify the adrenoleukodystrophy phenotype. It may also play a role in the process of peroxi


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ABCD4 Monoclonal Antibody