



ABCF2 Monoclonal Antibody

Catalog No	YP-mAb-00658
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABCF2
Immunogen	The antiserum was produced against synthesized peptide derived from human ABCF2. AA range:171-220
Specificity	ABCF2 Monoclonal Antibody detects endogenous levels of ABCF2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABC28;ATP-binding cassette sub-family F member 2;HUSY18;HUSY-18;Iron-inhibited ABC transporter 2;ABCF2
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	nucleosome,nucleus,mitochondrion,mitochondrial envelope,membrane,ATP-binding cassette (ABC) transporter complex,
Tissue Specificity	Colon,Lung,Uterus,
Function	caution:Lacks transmembrane domains and is probably not involved in transport.,similarity:Belongs to the ABC transporter family.,similarity:Belongs to the ABC transporter family. EF3 subfamily.,similarity:Contains 2 ABC transporter domains.,
Background	This gene encodes a member of the ATP-binding cassette (ABC) transporter superfamily. ATP-binding cassette proteins transport various molecules across extra- and intracellular membranes. Alterations in this gene may be involved in cancer progression. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 3 and 7. [provided



by RefSeq, Jul 2013],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ABCF2 Monoclonal Antibody