



ACBD6 Monoclonal Antibody

Catalog No	YP-mAb-00660
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACBD6
Immunogen	The antiserum was produced against synthesized peptide derived from human ACBD6. AA range:121-170
Specificity	ACBD6 Monoclonal Antibody detects endogenous levels of ACBD6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACBD6;Acyl-CoA-binding domain-containing protein 6
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Detected in placenta and spleen (at protein level). Detected in placenta, umbilical cord blood, CD34-positive hematopoietic progenitor cells and bone marrow.
Function	function:Binds long-chain acyl-coenzyme A molecules with a strong preference for unsaturated C18:1-CoA, lower affinity for unsaturated C20:4-CoA, and very weak affinity for saturated C16:0-CoA. Does not bind fatty acids., similarity:Contains 1 ACB (acyl-CoA-binding) domain.,similarity:Contains 2 ANK repeats.,subunit:Monomer.,tissue specificity:Detected in placenta and spleen (at protein level). Detected in placenta, umbilical cord blood, CD34-positive hematopoietic progenitor cells and bone marrow.,



Background

function: Binds long-chain acyl-coenzyme A molecules with a strong preference for unsaturated C18:1-CoA, lower affinity for unsaturated C20:4-CoA, and very weak affinity for saturated C16:0-CoA. Does not bind fatty acids., similarity: Contains 1 ACB (acyl-CoA-binding) domain., similarity: Contains 2 ANK repeats., subunit: Monomer., tissue specificity: Detected in placenta and spleen (at protein level). Detected in placenta, umbilical cord blood, CD34-positive hematopoietic progenitor cells and bone marrow.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ACBD6 Monoclonal Antibody