



AI-BP Monoclonal Antibody

Catalog No	YP-mAb-00663
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	APOA1BP
Immunogen	The antiserum was produced against synthesized peptide derived from human AIBP. AA range:121-170
Specificity	AI-BP Monoclonal Antibody detects endogenous levels of AI-BP protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AIBP;AI-BP;APOA1BP;Apolipoprotein A-I-binding protein;EC:5.1.99.6;NAD(P)H-hydrate epimerase;YJEFN1;NAD(P)HX epimerase;YjeF N-terminal domain-containing protein 1;YjeF_N1
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Mitochondrion . Secreted . In sperm, secretion gradually increases during capacitation. .
Tissue Specificity	Ubiquitously expressed, with highest levels in kidney, heart and liver. Present in cerebrospinal fluid and urine but not in serum from healthy patients. Present in serum of sepsis patients (at protein level).
Function	similarity:Contains 1 YjeF N-terminal domain.,subunit:Interacts with APOA1 and APOA2.,tissue specificity:Ubiquitously expressed, with highest levels in kidney, heart and liver. Present in cerebrospinal fluid and urine but not in serum from healthy patients. Present in serum of sepsis patients (at protein level).,
Background	The product of this gene interacts with apolipoprotein A-I (apoA-I), the major apolipoprotein of high-density lipoproteins (HDLs). It is secreted into some bodily fluids, and its synthesis and secretion are stimulated in vitro by incubating cells



with apoA-I. The human genome contains related pseudogenes. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using AI-BP Monoclonal Antibody