



ApoF Monoclonal Antibody


Catalog No	YP-mAb-00668
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	APOF
Immunogen	The antiserum was produced against synthesized peptide derived from human APOF. AA range:221-270
Specificity	ApoF Monoclonal Antibody detects endogenous levels of ApoF protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Apo F;Apo-F;apolipoprotein F;LTIP;ApoF;Lipid transfer inhibitor protein
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Secreted .
Tissue Specificity	Expressed by the liver and secreted in plasma.
Function	function:Minor apolipoprotein that associates with LDL. Inhibits cholesteryl ester transfer protein (CETP) activity and appears to be an important regulator of cholesterol transport. Also associates to a lesser degree with VLDL, Apo-AI and Apo-AII.,tissue specificity:Liver. Secreted in plasma.,
Background	The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ApoF Monoclonal Antibody