




ApoOL Monoclonal Antibody

Catalog No	YP-mAb-00674
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	APOOL
Immunogen	The antiserum was produced against synthesized peptide derived from human APOOL. AA range:219-268
Specificity	ApoOL Monoclonal Antibody detects endogenous levels of ApoOL protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	apolipoprotein O like;APOOL;CXorf33;FAM121A;MIC27;Protein FAM121A; UNQ8193;Apolipoprotein O-like
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein . Mitochondrion .
Tissue Specificity	Liver,
Function	similarity:Belongs to the apolipoprotein O family.,
Background	This gene encodes a protein which contains an apolipoprotein O superfamily domain. This domain is found on proteins in circulating lipoprotein complexes. [provided by RefSeq, Sep 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ApoOL Monoclonal Antibody