



StARD10 Monoclonal Antibody

Catalog No	YP-mAb-00736
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	STARD10
Immunogen	Synthesized peptide derived from the Internal region of human StARD10.
Specificity	StARD10 Monoclonal Antibody detects endogenous levels of StARD10 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Antigen NY CO 28;CGI 52;NY CO 28;PCTP L;PCTP like protein;PCTP2;SDCCAG28;STARD10;PCTP-like protein;CGI-52;PCTP-L;Antigen NY-CO-28;START domain-containing protein 10;Serologically defined colon cancer antigen 28;StAR-related lipid transfer protein 10
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	Cell projection, cilium, flagellum . Cytoplasm . Membrane . In testis was predominantly detected at the flagella of elongated spermatids, with a strong signal also found at the tail of epididymal sperm (By similarity). Mainly cytosolic. .
Tissue Specificity	Colon carcinoma,Uterus,
Function	function:May play specific roles in sperm maturation or fertilization.,sequence caution:Various sequencing problems as well as a translation in a wrong frame.,similarity:Contains 1 START domain.,
Background	function:May play specific roles in sperm maturation or fertilization.,sequence caution:Various sequencing problems as well as a translation in a wrong frame.,similarity:Contains 1 START domain.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using StARD10 Monoclonal Antibody