



ENT1 Monoclonal Antibody

Catalog No	YP-mAb-00757
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC29A1
Immunogen	The antiserum was produced against synthesized peptide derived from human ENT1. AA range:15-64
Specificity	ENT1 Monoclonal Antibody detects endogenous levels of ENT1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ENT-1;Equilibrative NBMPR-sensitive nucleoside transporter;Equilibrative;ENT1; SLC29A1;Equilibrative nucleoside transporter 1;Equilibrative nitrobenzylmercaptapurine riboside-sensitive nucleoside transporter;Nucleoside transporter;es-type;Solute carrier family 29 member 1; es-type
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Basolateral cell membrane; Multi-pass membrane protein. Apical cell membrane; Multi-pass membrane protein. Cell membrane; Multi-pass membrane protein. Predominantly localized in the basolateral membrane in polarized MDCK cells.
Tissue Specificity	Detected in erythrocytes (at protein level). Expressed in heart, brain, mammary gland, erythrocytes and placenta, and also in fetal liver and spleen.
Function	function:Mediates both influx and efflux of nucleosides across the membrane (equilibrative transporter). It is sensitive (ES) to low concentrations of the inhibitor nitrobenzylmercaptapurine riboside (NBMPR) and is sodium-independent. It has a higher affinity for adenosine. Inhibited by dipyrindamole and dilazep (anticancer chemotherapeutics drugs).,PTM:Glycosylated., similarity:Belongs to the SLC29A transporter family.,tissue specificity:Expressed in heart, brain, mammary gland, erythrocytes and placenta, and also in fetal liver



and spleen.,

Background

This gene is a member of the equilibrative nucleoside transporter family. The gene encodes a transmembrane glycoprotein that localizes to the plasma and mitochondrial membranes and mediates the cellular uptake of nucleosides from the surrounding medium. The protein is categorized as an equilibrative (as opposed to concentrative) transporter that is sensitive to inhibition by nitrobenzylthioinosine (NBMPR). Nucleoside transporters are required for nucleotide synthesis in cells that lack de novo nucleoside synthesis pathways, and are also necessary for the uptake of cytotoxic nucleosides used for cancer and viral chemotherapies. Multiple alternatively spliced variants, encoding the same protein, have been found for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ENT1 Monoclonal Antibody