



IGFN1 Monoclonal Antibody


Catalog No	YP-mAb-01204
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IGFN1
Immunogen	Recombinant Protein of IGFN1
Specificity	The antibody detects endogenous IGFN1 protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KYIP1;EEF1A2BP1;IGFN1;Immunoglobulin-like and fibronectin type III domain-containing protein 1;EEF1A2-binding protein 1;KY-interacting protein 1
Observed Band	115kD
Calculated Molecular Weight	115kD
Cell Pathway	Nucleus . Cytoplasm, myofibril, sarcomere, Z line .
Tissue Specificity	Expressed in skeletal muscle.
Function	similarity:Contains 4 fibronectin type-III domains.,similarity:Contains 5 Ig-like (immunoglobulin-like) domains.,subunit:Interacts with KY.,tissue specificity:Expressed in skeletal muscle.,
Background	similarity:Contains 4 fibronectin type-III domains.,similarity:Contains 5 Ig-like (immunoglobulin-like) domains.,subunit:Interacts with KY.,tissue specificity:Expressed in skeletal muscle.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using IGFBP1 Monoclonal Antibody