



eEF1A2 binding protein Monoclonal Antibody

Catalog No	YP-mAb-01212
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IGFN1
Immunogen	Recombinant Protein of eEF1A2 binding protein
Specificity	The antibody detects endogenous eEF1A2 binding protein protein
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	eEF1A2 binding protein;IGFN1;Immunoglobulin-like and fibronectin type III domain-containing protein 1 (EEF1A2-binding protein 1) (KY-interacting protein 1);Immunoglobulin-like and fibronectin type III domain-containing protein 1 (EEF1A2-binding protein 1;KY-interacting protein 1)
Observed Band	115kD
Calculated Molecular Weight	115kD
Cell Pathway	Nucleus . Cytoplasm, myofibril, sarcomere, Z line .
Tissue Specificity	Expressed in skeletal muscle.
Function	similarity:Contains 4 fibronectin type-III domains.,similarity:Contains 5 Ig-like (immunoglobulin-like) domains.,subunit:Interacts with KY.,tissue specificity:Expressed in skeletal muscle.,
Background	similarity:Contains 4 fibronectin type-III domains.,similarity:Contains 5 Ig-like (immunoglobulin-like) domains.,subunit:Interacts with KY.,tissue specificity:Expressed in skeletal muscle.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

130KD —
95KD —
70KD —
55KD —



Western Blot analysis of various cells using eEF1A2 binding protein Monoclonal Antibody