



CACNG2 Monoclonal Antibody

Catalog No	YP-mAb-01217
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CACNG2
Immunogen	Synthetic Peptide of CACNG2 AA range: 58-108
Specificity	The antibody detects endogenous CACNG2 protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CACNG2;Voltage-dependent calcium channel gamma-2 subunit;Neuronal voltage-gated calcium channel gamma-2 subunit;Transmembrane AMPAR regulatory protein gamma-2;TARP gamma-2
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Membrane; Multi-pass membrane protein. Cell junction, synapse, synaptosome .
Tissue Specificity	Brain.
Function	function:Thought to stabilize the calcium channel in an inactivated (closed) state., PTM:Phosphorylation of Thr-321 impairs interaction with DLG1 and DLG4., similarity:Belongs to the PMP-22/EMP/MP20 family. CACNG subfamily., subunit:The L-type calcium channel is composed of five subunits: alpha-1, alpha-2/delta, beta and gamma. Interacts with the PDZ domains of DLG4/PSD-95 and DLG1/SAP97. May interact with GOPC.,tissue specificity:Brain.,
Background	calcium voltage-gated channel auxiliary subunit gamma 2(CACNG2) Homo sapiens The protein encoded by this gene is a type I transmembrane AMPA receptor regulatory protein (TARP). TARPs regulate both trafficking and channel



gating of the AMPA receptors. This gene is part of a functionally diverse eight-member protein subfamily of the PMP-22/EMP/MP20 family. This gene is a susceptibility locus for schizophrenia. [provided by RefSeq, Dec 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CACNG2 Monoclonal Antibody