




Brp16 Monoclonal Antibody

Catalog No	YP-mAb-01572
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM203A
Immunogen	The antiserum was produced against synthesized peptide derived from human BRP16. AA range:251-300
Specificity	Brp16 Monoclonal Antibody detects endogenous levels of Brp16 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C8orf30A;HGH1;Brain protein 16;Brp16;FAM203A;Protein FAM203A
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	
Tissue Specificity	Colon,Lung,Umbilical cord blood,
Function	similarity:Belongs to the UPF0570 (BRP16/HGH1) family.,
Background	similarity:Belongs to the UPF0570 (BRP16/HGH1) family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Brp16 Monoclonal Antibody