




BTBD6 Monoclonal Antibody

Catalog No	YP-mAb-01576
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BTBD6
Immunogen	The antiserum was produced against synthesized peptide derived from human BTBD6. AA range:91-140
Specificity	BTBD6 Monoclonal Antibody detects endogenous levels of BTBD6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BDPL;BTB (POZ) domain containing 6;BTBD6;Lens BTB domain protein;BTB/POZ domain-containing protein 6
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Cytoplasm . Found in punctated bodies in the cytoplasm. .
Tissue Specificity	Expressed in lens.
Function	similarity:Contains 1 BTB (POZ) domain.,tissue specificity:Expressed in lens.,
Background	similarity:Contains 1 BTB (POZ) domain.,tissue specificity:Expressed in lens.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using BTBD6 Monoclonal Antibody