




BUD31 Monoclonal Antibody

Catalog No	YP-mAb-01578
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BUD31
Immunogen	The antiserum was produced against synthesized peptide derived from human BUD31. AA range:10-59
Specificity	BUD31 Monoclonal Antibody detects endogenous levels of BUD31 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BUD31;BUD31 homolog (S. cerevisiae);EDG 2;EDG2;fSAP17;G10;Protein BUD31 homolog;Protein EDG 2;Protein G10 homolog;YCR063W;Protein EDG-2
Observed Band	17kD
Calculated Molecular Weight	17kD
Cell Pathway	Nucleus . Detected in chromatin at the promoter of AR target genes. .
Tissue Specificity	Detected in epithelial and stromal cells in benign prostate hyperplasia tissue (at protein level).
Function	similarity:Belongs to the BUD31 (G10) family.,
Background	similarity:Belongs to the BUD31 (G10) family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using BUD31 Monoclonal Antibody