




CCDC4 Monoclonal Antibody

Catalog No	YP-mAb-01597
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BEND4
Immunogen	The antiserum was produced against synthesized peptide derived from human BEND4. AA range:371-420
Specificity	CCDC4 Monoclonal Antibody detects endogenous levels of CCDC4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BEN domain containing 4;BEND4;CCDC4;BEN domain-containing protein 4; Coiled-coil domain-containing protein 4
Observed Band	57kD
Calculated Molecular Weight	57kD
Cell Pathway	
Tissue Specificity	Prostatic carcinoma,Testis,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 BEN domain.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 BEN domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CCDC4 Monoclonal Antibody