




CCDC45 Monoclonal Antibody

Catalog No	YP-mAb-01598
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CEP95
Immunogen	The antiserum was produced against synthesized peptide derived from human CCDC45. AA range:591-640
Specificity	CCDC45 Monoclonal Antibody detects endogenous levels of CCDC45 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCDC45;CEP95;DKFZp667E1824;Centrosomal protein of 95 kDa;CEP45;Coiled-coil domain-containing protein 45
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole .
Tissue Specificity	Lymph,Lymph node,Testis,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CCDC45 Monoclonal Antibody