




CGI-43 Monoclonal Antibody

Catalog No	YP-mAb-01608
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCZ1/C7orf28B
Immunogen	The antiserum was produced against synthesized peptide derived from human CG028. AA range:411-460
Specificity	CGI-43 Monoclonal Antibody detects endogenous levels of CGI-43 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C7orf28A;C7orf28B;CGI 43;Chromosome 7 open reading frame 28A;Chromosome 7 open reading frame 28B;DKFZp586I1023;FLJ60592;H_DJ1163J12.2;H_NH0577018.2;Hypothetical protein LOC221960;Hypothetical protein LOC51622;MGC19819;UPF0550 protein C7orf28;CGI-43;CCZ1/C7orf28B;Vacuolar fusion protein CCZ1 homolog/Vacuolar fusion protein CCZ1 homolog-like;CCZ1;Vacuolar fusion protein CCZ1 homolog;CCZ1B;Vacuolar fusion protein CCZ1 homolog B;Vacuolar fusion protein CCZ1 homolog-like
Observed Band	56kD
Calculated Molecular Weight	56kD
Cell Pathway	Lysosome membrane .
Tissue Specificity	Brain,Skin,Testis,Uterus,
Function	similarity:Belongs to the UPF0550 family.,
Background	similarity:Belongs to the UPF0550 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CGI-43 Monoclonal Antibody