



CRP1 Monoclonal Antibody


Catalog No	YP-mAb-01623
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CSRP1
Immunogen	The antiserum was produced against synthesized peptide derived from human CRP1. AA range:51-100
Specificity	CRP1 Monoclonal Antibody detects endogenous levels of CRP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRP-1;CRP1;CRP 1;CRIP 1;CRIP;CSRP1;Cysteine and glycine-rich protein 1; CSRP;CYRP;Cysteine-rich protein 1;CRP
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Nucleus .
Tissue Specificity	Amygdala,Brain,Cervix,Epithelium,Pancreas,Term placenta,Trachea,Ute
Function	function:Could play a role in neuronal development.,similarity:Contains 1 LIM zinc-binding domain.,similarity:Contains 2 LIM zinc-binding domains.,subunit:Monomer.,
Background	This gene encodes a member of the cysteine-rich protein (CSRP) family. This gene family includes a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. The LIM/double zinc-finger motif found in this gene product occurs in proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2010],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CRP1 Monoclonal Antibody