



CRSP130 Monoclonal Antibody

Catalog No	YP-mAb-01624
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MED23
Immunogen	The antiserum was produced against synthesized peptide derived from human MED23. AA range:1-50
Specificity	CRSP130 Monoclonal Antibody detects endogenous levels of CRSP130 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MED23/CRSP3;CRSP133;CRSP130;CRSP complex subunit 3;Cofactor required for Sp1 transcriptional activation subunit 3;MED23;Mediator of RNA polymerase II transcription subunit 23;ARC130;CRSP3;DRIP130;KIAA1216;SUR2;Activator-recruited cofactor 130 kDa component;CRSP complex subu
Observed Band	150kD
Calculated Molecular Weight	150kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Cervix carcinoma,Placenta,Prostate,Testis,
Function	alternative products:Experimental confirmation may be lacking for some isoforms,function:Required for transcriptional activation subsequent to the assembly of the preinitiation complex (By similarity). Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general



transcription factors. Required for transcriptional activation by adenovirus E1A protein. Required for ELK1-dependent transcriptional activation in response to activated Ras signaling. similarity: Belongs to the

Background

The activation of gene transcription is a multistep process that is triggered by factors that recognize transcriptional enhancer sites in DNA. These factors work with co-activators to direct transcriptional initiation by the RNA polymerase II apparatus. The protein encoded by this gene is a subunit of the CRSP (cofactor required for SP1 activation) complex, which, along with TFIID, is required for efficient activation by SP1. This protein is also a component of other multisubunit complexes e.g. thyroid hormone receptor-(TR-) associated proteins which interact with TR and facilitate TR function on DNA templates in conjunction with initiation factors and cofactors. This protein also acts as a metastasis suppressor. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CRSP130 Monoclonal Antibody